


INSTALLATION INSTRUCTIONS

The Battery Back-up System allows the Guardian Alert 911 product to continue to function in the event of a power failure. This product is an accessory to the Guardian Alert 911 product.

The system comes with 6 NiMH AA batteries (1500mAh).

Installation:

1. If the batteries are wrapped in plastic wrap - please unwrap all 6 batteries.
2. Install the 6 batteries that came with the product into the Battery Back-up module. To do this, slide the blue cover forward until the battery cover comes off. Install the 6 batteries with the positive terminals for all 6 batteries facing the front. Slide the blue cover back on to the unit.
3. Unplug the 12 volt power adaptor that normally goes to the back of the Guardian Alert 911 base unit - and put the Battery Back-up System in between. In other words - plug the 12 volt power adaptor into the Battery Back-up System - and plug the cord from the Battery Back-Up System into the back of the base station. This is shown in the diagram above.
4. Plug the 12 volt power adaptor into the wall power outlet. Be sure to pick a power outlet that has power continually - one that is not turned off with a light switch. At this time, the green LED light should come on - and stay on for as long as you have the power adaptor plugged into the wall outlet.



The red LED light on the unit will come on during charging. Once the batteries are charged, the red LED light will go out.

NOTE: It is recommended that you change the batteries every two to three years. When getting new batteries - use NiMH AA batteries of 1,500 mAh or higher. With 1500 mAh batteries - the unit should power the Guardian Alert 911 system for 24 hours. Higher amp-hour batteries (i.e. 2500 mAh) will power the units almost twice as long.

CAUTION: Periodically check that your 12 volt wall adaptor is firmly plugged into the wall outlet. The fact that the green LED light is on signifies this to be the case.

LogicMark

8625 Hampton Way; Fairfax Station, Virginia 22039

Tel: 703-934-7934

www.LogicMark.com